

Product	ProMix™ Survivin Peptide Pool
No. of Peptides	11 peptides
Quantity	approx. 25 µg / peptide
Description	This pool consists of 11 peptides, each corresponding to a defined HLA class I-restricted T-cell epitope from Survivin.
Sequences	<i>see overleaf</i>
Purity (HPLC)	>90% (determined at 220 nm; individual peptides)
Remark	The peptides of this product are supplied as trifluoroacetate salts
Storage	Store at - 20°C
Solubilization	Dissolve in a minimum amount of pure DMSO (approx. 40 µl) and dilute with PBS buffer to the final concentration. Please note that the final concentration of DMSO must be below 1 % (v/v) to avoid toxicity in a biological system.

Use recommended within 6 months.

Released by:
(dated as per product label above)

Sequences:

Peptide	Sequence	Epitope Source	HLA Allele Restriction
1	MAEAGFIHY	Survivin ₃₈₋₄₆ (mod)	A*01:01
2	PTENEPDLAY	Survivin ₄₇₋₅₆ (mod)	A*01:01
3	ELTLGEFLKL	Survivin ₉₅₋₁₀₄	A*02:01
4	LMLGEFLKL	Survivin ₉₆₋₁₀₄ (mod)	A*02:01
5	TLPPAWQPFL	Survivin ₅₋₁₄	A*02:01
6	RISTFKNWP	Survivin ₁₈₋₂₇ (mod)	A*03:01
7	DLAQCFEFCFK	Survivin ₅₆₋₆₂	A*11:01
8	AYACNTSTL	Survivin ₈₀₋₈₈	A*24:02
9	STFKNWPFL	Survivin ₂₀₋₂₈	A*24:02
10	LPPAWQPFL	Survivin ₆₋₁₄	B*07:02
11	EPDLAQCFY	Survivin ₅₁₋₅₉	B*35:01

© Copyright ProImmune Limited 2007-2022. All Rights Reserved